

WORKING ELEMENT

passive cutting; closed handle

for Ø 2.9mm telescopes; HF-connection up

for mono- and bipolar use; ATTENTION: our bipolar resectoscopy system can be used only with HF-units which have the possibility of bipolar cut in saline solution!!



SKU: 630-335-60B

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 2.9
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 2.9
TYPE	Working Elements