

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

**040-580-24**


## **ROGERS Hysterectomy forceps**

straight; serrated

25.5cm (10")

Categories: Open Surgery, Hemostats

1 pcs.



## **ADDITIONAL INFORMATION**

<b>UNIT</b>	pcs.
<b>TYPE</b>	Hysterectomy Forceps
<b>Aesculap</b>	BJ550
<b>KLS Martin</b>	13-637-10
<b>AUTHORS</b>	ROGERS
<b>Medicon</b>	16.35.50
<b>SPECIFICATIONS</b>	serrated, straight