

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

040-581-24


ROGERS Hysterectomy forceps

slightly curved; serrated

25.0cm (10")

Categories: Open Surgery, Hemostats

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ551
KLS Martin	13-637-20
AUTHORS	ROGERS
Medicon	16.35.51
SPECIFICATIONS	curved, serrated