

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

040-583-21


ROGERS Hysterectomy forceps

curved; serrated

21.0cm (8.6")

Categories: Open Surgery, Hemostats

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ532
KLS Martin	13-635-30
AUTHORS	ROGERS
SPECIFICATIONS	curved, serrated