

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

**040-583-24**


## **ROGERS Hysterectomy forceps**

curved; serrated

24.5cm (9.6")

Categories: Open Surgery, Hemostats

1 pcs.



## **ADDITIONAL INFORMATION**

<b>UNIT</b>	pcs.
<b>TYPE</b>	Hysterectomy Forceps
<b>Aesculap</b>	BJ552
<b>KLS Martin</b>	13-637-30
<b>AUTHORS</b>	ROGERS
<b>SPECIFICATIONS</b>	curved, serrated