

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

040-585-21


ROGERS Hysterectomy forceps

angled; serrated

21.0cm (8 1/2")

Categories: Open Surgery, Hemostats

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ533
KLS Martin	13-635-40
AUTHORS	ROGERS
SPECIFICATIONS	angled, serrated