

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

**040-585-24**


## **ROGERS Hysterectomy forceps**

angled; serrated

24.0cm (9.6")

Categories: Open Surgery, Hemostats

1 pcs.



## **ADDITIONAL INFORMATION**

<b>UNIT</b>	pcs.
<b>TYPE</b>	Hysterectomy Forceps
<b>Aesculap</b>	BJ553
<b>KLS Martin</b>	13-637-40
<b>AUTHORS</b>	ROGERS
<b>Medicon</b>	16.35.53
<b>SPECIFICATIONS</b>	angled, serrated