

TONTARRA / Open Surgery / Hemostats / Hysterectomy Forceps / ROGERS Hysterectomy forceps

040-580-21


ROGERS Hysterectomy forceps

straight; serrated

21.5cm (8.6")

Categories: Open Surgery, Hemostats

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ530
KLS Martin	13-635-10
AUTHORS	ROGERS
SPECIFICATIONS	serrated, straight