

TONTARRA / Open Surgery / Retractors / Selfretaining Retractors / FINSEN selfretaining Retractor

070-610-07


FINSEN selfretaining Retractor

sharp; 3X4 teeth

7.0cm (2.8")

Categories: Open Surgery, Retractors

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV002
KLS Martin	15-702-07
AUTHORS	FINSEN
Medicon	24.04.07
SPECIFICATIONS	sharp, toothed