

TONTARRA / Open Surgery / Retractors / Selfretaining Retractors / FINSEN selfretaining Retractor

**070-611-05**


## **FINSEN selfretaining Retractor**

blunt; 2X3 teeth

5.0cm (2")

Categories: Open Surgery, Retractors

1 pcs.



### **ADDITIONAL INFORMATION**

|                       |                          |
|-----------------------|--------------------------|
| <b>UNIT</b>           | pcs.                     |
| <b>TYPE</b>           | Selfretaining Retractors |
| <b>KLS Martin</b>     | 15-701-05                |
| <b>AUTHORS</b>        | FINSEN                   |
| <b>Medicon</b>        | 24.04.06                 |
| <b>SPECIFICATIONS</b> | blunt, toothed           |