

TONTARRA / Open Surgery / Retractors / Selfretaining Retractors / FINSEN selfretaining Retractor

070-611-07


FINSEN selfretaining Retractor

blunt; 3X4 teeth

7.0cm (2.8")

Categories: Open Surgery, Retractors

1 pcs.



ADDITIONAL INFORMATION

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV003
KLS Martin	15-703-07
AUTHORS	FINSEN
Medicon	24.04.08
SPECIFICATIONS	blunt, toothed