

2- LUMEN RETRIEVAL BASKET

hexagonal, 30 x 60mm; 4 wires

3.3mm; 210cm; for guide wire 0.035"



SKU: PK-1FHPK3330230K-2L

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)