

## 3-RING SYRINGE WITH GLASS BARREL

Luer-Lock connection; 10 ml



**SKU:** 087-055-10

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Syringes</a>
<b>KLS Martin</b>	18-915-10
<b>Medicon</b>	09.78.54

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Syringes</a>
<b>KLS Martin</b>	18-915-10
<b>Medicon</b>	09.78.54