

## 3-WAY STOPCOCK

1X LL= male; 2X LL= female



**SKU:** 087-922-00

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stopcocks and Connecting Pieces</a>
<b>Aesculap</b>	SK320
<b>KLS Martin</b>	18-845-00
<b>Medicon</b>	09.79.65

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stopcocks and Connecting Pieces</a>
<b>Aesculap</b>	SK320
<b>KLS Martin</b>	18-845-00
<b>Medicon</b>	09.79.65