

## 45° ANGLED NEPHROSCOPE

19Fr., working channel Ø 5.0mm

view 15°, working length 200mm

**SKU:** 601-130-45

**Categories:** [Endoscopy](#), [Equipment](#), [Urology](#)




### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Urology](#)