


ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth

16.5cm (6.6")



SKU: 070-860-16

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV995
KLS Martin	15-771-16
AUTHORS	ADSON
Medicon	24.30.16
SPECIFICATIONS	sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV995
KLS Martin	15-771-16

AUTHORS

[ADSON](#)

Medicon

24.30.16

SPECIFICATIONS

[sharp](#), [toothed](#)