


ADSON (CONE) SELFRETAINING RETRACTOR

sharp; 3X4 teeth
 20.0cm (8")



SKU: 070-860-20
Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-770-20
AUTHORS	ADSON
Medicon	24.30.20
SPECIFICATIONS	sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-770-20
AUTHORS	ADSON
Medicon	24.30.20

SPECIFICATIONS

[sharp](#), [toothed](#)