


ADSON (CONE) SELFRETAINING RETRACTOR

blunt; 3X4 teeth

26.0cm (10.4")



SKU: 070-861-26
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-771-26
AUTHORS	ADSON
SPECIFICATIONS	blunt , toothed