

**ADSON DRESSING FORCEPS**

straight; cross serrated


12.0cm (4.8"); width of jaw 1.3 +/-0.1

**SKU:** 020-078-12  
**Categories:** [Open Surgery](#), [Forceps](#)



**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                                   |
|----------------|-------------------------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                                              |
| TYPE           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| AUTHORS        | <a href="#">CUSHING</a>                                           |
| Medicon        | 06.20.14                                                          |
| SPECIFICATIONS | <a href="#">serrated</a> , <a href="#">straight</a>               |



**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                                   |
|----------------|-------------------------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                                              |
| TYPE           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| AUTHORS        | <a href="#">CUSHING</a>                                           |
| Medicon        | 06.20.14                                                          |
| SPECIFICATIONS | <a href="#">serrated</a> , <a href="#">straight</a>               |