

ADSON NEEDLE HOLDER

straight; serrated; fenestrated

18.0cm (7.2")




SKU: 050-192-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM224
KLS Martin	20-360-18
AUTHORS	ADSON
Medicon	10.36.17
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM224
KLS Martin	20-360-18

AUTHORS

[ADSON](#)

Medicon

10.36.17

SPECIFICATIONS

[serrated](#), [straight](#)