

## ADSON NEEDLE HOLDER

straight; serrated; fenestrated

18.0cm (7.2")




**SKU:** 050-192-18

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM224
<b>KLS Martin</b>	20-360-18
<b>AUTHORS</b>	<a href="#">ADSON</a>
<b>Medicon</b>	10.36.17
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM224
<b>KLS Martin</b>	20-360-18

## AUTHORS

[ADSON](#)

## Medicon

10.36.17

## SPECIFICATIONS

[serrated](#), [straight](#)