

ADSON SUCTION TUBE

11 Fr.; angled

14.5cm (5.8") working length

SKU: 085-050-11

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | ADSON |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | ADSON |