

## **ALLIGATOR FORCEPS; FL 600MM**

5 Fr. semi-rigid

double action; insert only



**SKU:** 583-460-16A  
**Categories:** [Endoscopy](#), [Urology](#)

### **ZUSÄTZLICHE INFORMATIONEN**

|       |                                        |
|-------|----------------------------------------|
| UNIT  | <a href="#">pcs.</a>                   |
| RANGE | <a href="#">Uretero-Renoscscopy</a>    |
| TYPE  | <a href="#">Semi-rigid Instruments</a> |



## **ZUSÄTZLICHE INFORMATIONEN**

|       |                                        |
|-------|----------------------------------------|
| UNIT  | <a href="#">pcs.</a>                   |
| RANGE | <a href="#">Uretero-Renoscscopy</a>    |
| TYPE  | <a href="#">Semi-rigid Instruments</a> |