


ALM SELFRETAINING RETRACTOR

semi-sharp; 4X4 teeth

10.0cm (4")



SKU: 070-622-10

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ALM
SPECIFICATIONS	semi-sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ALM
SPECIFICATIONS	semi-sharp , toothed