

ANDERSON-ADSON SELFRETAINING RETRACTOR

blunt; 4X4 teeth; curved


19.0cm (7.6")



SKU: 070-679-20
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-732-20
AUTHORS	ADSON , ANDERSON
Medicon	24.28.02
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-732-20
AUTHORS	ADSON , ANDERSON
Medicon	24.28.02

SPECIFICATIONS

[sharp](#), [toothed](#)