

ANDREWS-PYNCHON SUCTION TUBE

removable suction tip

inner Ø 2.0mm; outer Ø 3.0mm; 24.0cm

SKU: 085-006-27

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | ANDREWS-PYNCHON |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | ANDREWS-PYNCHON |