

ARTHROSCOPE/ SINUSCOPE HD+

Ø 4.0mm; 70°; 175mm

autoclavable




SKU: 601-204-70HD

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	Arthroscopy , ENT



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	Arthroscopy , ENT