

ARTHROSCOPE/ SINUSCOPE HD+

Ø 2.7mm; 70°; 140 mm

autoclavable

SKU: 601-201-70HD

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)