

## ARTHROSCOPE/ SINUSCOPE HD+

Ø 4.0mm; 45°; 175mm

autoclavable

**SKU:** 601-204-40HD

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Arthroscopy](#), [ENT](#)



45°

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Arthroscopy](#), [ENT](#)