

ARTHROSCOPE/ SINUSCOPE HD+2

Ø 4.0mm; 45°; 175mm

autoclavable

SKU: 601-204-40HD2

Categories: [Endoscopy](#), [Equipment](#), [Arthroscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)