

ARTHROSCOPE/ SINUSCOPE HD+2

Ø 4.0mm; 70°; 175mm

autoclavable

SKU: 601-204-70HD2

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)