

ARTHROSCOPY TELESCOPE HD

Ø 2.7mm; 30°; 110 mm

autoclavable

SKU: 601-303-30

Categories: [Endoscopy](#), [Equipment](#), [Arthroscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)