

ARTHROSCOPY TELESCOPE HD

Ø 2.7mm; 0°; 110 mm

autoclavable

SKU: 601-303-00

Categories: [Endoscopy](#), [Equipment](#), [Arthroscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Arthroscopy](#), [ENT](#)