

ARTHROSCOPY TELESCOPE HD

Ø 2.7mm; 0°; 110 mm

autoclavable

SKU: 601-303-00

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Arthroscopy](#), [ENT](#)