

AUTOMATIC TROCAR -SLEEVE, Ø 11MM, WL 75MM

short, Stainless Steel, with Stopcock

and rough tube, distal bevelling



SKU: BLUO/04-00413-2-SF

Categories: [Endoscopy](#), [Blue endo](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)