

AUTOMATIC TROCAR -SLEEVE, Ø 13 MM,

WL120mm, stainless steel, maintenance free

stopcock, rough tube, distal bevelling



SKU: BLUO/04-00411-2-F

Categories: [Blue endo](#), [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Light Solution](#)