

## AUTOMATIC TROCAR -SLEEVE, Ø 5.5MM, WL75MM

short, Stainless Steel, maintenance-free

stopcock, rough tube, distal bevelling



**SKU:** BLUO/04-00415-2-SF

**Categories:** [Blue endo](#), [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Light Solution](#)