

BABY-ADSON HAEMOSTATIC FORCEPS

strongly curved; serrated

18.0cm (7.2")




SKU: 040-411-18

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
Aesculap	BJ014
KLS Martin	13-421-18
AUTHORS	ADSON , BABY
Medicon	17.05.19
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
Aesculap	BJ014
KLS Martin	13-421-18

AUTHORS

[ADSON](#), [BABY](#)

Medicon

17.05.19

SPECIFICATIONS

[curved](#), [serrated](#)