

BABY-DOYEN INTESTINAL CLAMP

straight; diagonally serrated

18.0cm (7.2")



SKU: 160-200-18

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA120
KLS Martin	30-276-18
AUTHORS	BABY , DOYEN
Medicon	50.16.17
SPECIFICATIONS	diagonally serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA120
KLS Martin	30-276-18

AUTHORS

[BABY](#), [DOYEN](#)

Medicon

50.16.17

SPECIFICATIONS

[diagonally serrated](#), [straight](#)