

BABY-DOYEN INTESTINAL CLAMP

straight; diagonally serrated

16.5cm (6.6")



SKU: 160-200-16

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA120
AUTHORS	BABY , DOYEN
SPECIFICATIONS	diagonally serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA120
AUTHORS	BABY , DOYEN
SPECIFICATIONS	diagonally serrated , straight