

BABY-MIKULICZ PERITONEAL CLAMP

curved; serrated; 1X2 teeth

14.0cm (5.6")




SKU: 040-701-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Peritoneal Clamp Forceps
Aesculap	BJ300
KLS Martin	13-721-14
AUTHORS	BABY , MIKULICZ
Medicon	16.85.14
SPECIFICATIONS	curved , serrated , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Peritoneal Clamp Forceps
Aesculap	BJ300
KLS Martin	13-721-14

AUTHORS

[BABY](#), [MIKULICZ](#)

Medicon

16.85.14

SPECIFICATIONS

[curved](#), [serrated](#), [toothed](#)