

## BAJONETT DISC RONGEUR

angled up; 4 x 10mm; 30mm

16.0cm (6.4")



**SKU:** 240-591-04

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">BAJONETT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">BAJONETT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a>