

BAJONETT DISC RONGEUR

angled up; 4 x 10mm; 30mm

16.0cm (6.4"); CC-coated



SKU: 240-591-04CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	angled up



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	angled up