

BAJONETT PUNCH; 40° UP CUT

5mm; without ejector

16.0cm (6.4"); 10mm open; CC-coated




SKU: 240-090-05CC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut, without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut, without ejector