

## BAJONETT PUNCH; 40° UP CUT

3mm; without ejector; thin footplate

16.0cm (6.4"); 10mm open




**SKU:** 240-090-03T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>AUTHORS</b>        | <a href="#">BAJONETT</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>   |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                          |
| <b>AUTHORS</b>        | <a href="#">BAJONETT</a>                                     |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">without ejector</a> |