

## BALLENGER SWIVEL KNIVES

bayonet; 3.0mm

21.0cm (8.2")

**SKU:** 112-231-03

**Categories:** [Open Surgery](#), [Rhinology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Septum Knives</a>
<b>Aesculap</b>	OL354
<b>KLS Martin</b>	37-505-03
<b>AUTHORS</b>	<a href="#">BALLENGER</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Septum Knives</a>
<b>Aesculap</b>	OL354
<b>KLS Martin</b>	37-505-03
<b>AUTHORS</b>	<a href="#">BALLENGER</a>