

BALLON-CATHETER

able for guide wire 0.035", Ø ballon 9mm

2.3mm, length 250cm, single use



SKU: PK-3LRPK5407250-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)