

## BARON SUCTION TUBE

3 Fr.; WL: 75mm



**SKU:** 086-300-10

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-528-10
<b>AUTHORS</b>	<a href="#">BARON</a>
<b>Medicon</b>	62.22.10

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-528-10
<b>AUTHORS</b>	<a href="#">BARON</a>
<b>Medicon</b>	62.22.10