

BECKMANN-ADSON SELFRETAINING RETRACTOR

blunt; 4X4 teeth; 25X31mm


31.0cm (12.4")



SKU: 070-871-02
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON , BECKMANN
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON , BECKMANN
SPECIFICATIONS	blunt , toothed