

BECKMANN-ADSON SELFRETAINING RETRACTOR

blunt; 4X4 teeth; 25X45mm


31.0cm (12.4")



SKU: 070-871-04
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-778-02
AUTHORS	ADSON , BECKMANN
Medicon	24.34.46
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-778-02
AUTHORS	ADSON , BECKMANN
Medicon	24.34.46

SPECIFICATIONS

[blunt](#), [toothed](#)